


Applying Forceps

As important as the aneurysm clips, the applying forceps has primordial part in the surgical act. On confirming this information, the standard ASTM F 700-93, at the article 5.1, advice: "Open and close an aneurysm clip with the clip applier specifically manufactured to handle that particularly aneurysm clip".

Vicca presents its applying forceps line for intracranial aneurysm clips.

Applying Forceps for Intracranial Aneurysm Clips


Bayonet



Availability:




Straight



Availability:




Short Bayonet



Availability:





Remover



Availability:





AVAILABILITY TABLE OF VICCA APPLYING FORCEPS

APPLYING FORCEPS FOR STANDARD CLIPS

VICCA CODE	TYPE			LENGTH			LOCK		TIP		
	BAYONET	BAYONET ARTICULATE	STRAIGHT	19 cm	22 cm	26 cm	WITH LOCK	WITHOUT LOCK	STRAIGHT	15°	MULTIANGLE
SG	✓			✓				✓	✓		
SH	✓			✓				✓		✓	
SD			✓	✓				✓	✓		
SB			✓	✓				✓		✓	
SL	✓				✓		✓	✓	✓		
SR	✓				✓		✓	✓		✓	
VS	✓				✓		✓	✓			✓
SE			✓		✓			✓	✓		
SF			✓		✓			✓		✓	
PS		✓			✓		✓	✓			✓
SO	✓					✓	✓	✓	✓		
SQ	✓					✓	✓	✓	✓	✓	

APPLYING FORCEPS FOR MINI CLIPS

VICCA CODE	TYPE			LENGTH			LOCK		TIP		
	BAYONET	BAYONET ARTICULATE	STRAIGHT	19 cm	22 cm	26 cm	WITH LOCK	WITHOUT LOCK	STRAIGHT	15°	MULTIANGLE
MG	✓			✓				✓	✓		
MH	✓			✓				✓		✓	
MD			✓	✓				✓	✓		
MB			✓	✓				✓		✓	
ML	✓				✓		✓	✓	✓		
MR	✓				✓		✓	✓		✓	
VM	✓				✓		✓	✓			✓
ME			✓		✓			✓	✓		
MF			✓		✓			✓		✓	
PM		✓			✓		✓	✓			✓
MO	✓					✓	✓	✓	✓		
MQ	✓					✓	✓	✓	✓	✓	

REMOVERS

VICCA CODE	PATTERN		TYPE			LENGTH			LOCK		TIP		
	STANDARD	MINI	BAYONET	BAYONET ARTICULATE	STRAIGHT	19 cm	22 cm	26 cm	WITH LOCK	WITHOUT LOCK	STRAIGHT	15°	MULTIANGLE
ES	✓		✓				✓			✓	✓		
EM		✓	✓				✓			✓	✓		

All applying forceps listed above are available for aneurysm clips in titanium and cobalt alloys.



Certified company

